SHERAprint-tray clear

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INSTRUCTIONS FOR USE

Description:

SHERAprint-tray clear is a light-curing resin for the additive fabrication of individual impression trays. The material has been optimized for Rapidshape DLP® systems with a wavelength of 385 nm for layer thicknesses of 150 µm.

Composition:

Formulation of acrylates and initiators. Please always observe the current safety data sheet!

Indication:

Additive fabrication of individual impression trays, bite blocks and set-up bases.

Contraindications:

Patient and user: products made from SHERAprint-tray clear should not be used if an allergy to one of the constituents is known. Otherwise no restrictions.

Adverse effects:

Patient and user: product may cause allergic reactions.

User: observe warnings. Unpolymerized material can cause skin irritation and severe eye irritation.

Interactions:

Interactions are not known.

Warnings:

Avoid inhaling vapors. Avoid contact with the uncured material and wear protective gloves, protective clothing and eye protection. In case of contact with the skin, wash off with copious amounts of water and soap. Seek medical advice in case of skin irritation or rashes.

Unpolymerized material is harmful to the environment; avoid releasing into the environment.

Please always observe the applicable country-specific occupational health and safety regulations as well as the legal requirements when handling hazardous substances.

Notes on storage / shelf life:

Protect from light. Can polymerize spontaneously. Store between 4°C [39°F] and 25°C [77°F]. Do not use after expiration date. See label for shelf life.

Notes on disposal:

Dispose of contents/ container in accordance with local regulations.

Application:

1. Design process

- Preparing the data set: the type of support structures can be found in the instructions for use of the system supplier. Possible alignment of the supporting structures:

labial or

occlusal.

- Defining variable process parameters: the selection of suitable process parameters is made according to the specifications of the system supplier. Care must be taken to select the correct material-specific parameter set (product, layer thickness, exposure time)

- Pouring the material: stir the material thoroughly with a spatula and pour into the material tray of the device. Then wait for 10 minutes.

- Start the design process.

2. Clean and dry

- Option 1: Cleaning the assembly: cleaning is performed under flowing conditions with an ultrasonic cleaner for 3D printed resin objects (e.g. SHERAultra-p) for 2,5 minutes. If necessary, clean the objects carefully with a brush and a spray bottle filled with an ultrasonic cleaner.

- Cleaning the assembly for a second time: we recommend cleaning the assembly for a second time in fresh ultrasonic cleaner following the above mentioned steps.

- Note: the maximum cleaning time for both steps combined is 5 minutes.

- Option 2: automatic cleaning in the SHERAprint-wash. Follow the instructions for the device.

- Drying the assembly: blow off assembly with compressed air and allow to dry for 30 minutes. As an alternative, the objects can be stored in a drying cabinet at a temperature of approx. 37 °C.

- If necessary, remove support structures

3. Post-exposure

- For post-exposure, the assemblies must be dry and have cooled to room temperature.

- Option 1: Post-exposure in the SHERAflash-light plus (Otoflash G17L, [NK-Optik]) or an identical unit with 2000 flashes each for the top and bottom surfaces of the object (Σ4000).

- Recommendation: the inhibition layer is prevented during curing under exclusion of oxygen.

- Option 2: Post-exposure in SHERAprint-cure with selection of the corresponding parameter set of the material. [Exposure time: 3 minutes / vacuum 50 mbar / without nitrogen / heating off / power 75 % / wavelength 375+415 nm].

- Care must be exercised to maintain or calibrate the equipment used in accordance with the manufacturer's specifications to guarantee the mechanical properties, the color and the biocompatibility.

4. Additional notes on processing:

Assemblies made from SHERAprint-tray clear are subsequently processed and finished according to the usual dental technical aspects.

Unit of delivery:

SHERAprint-tray clear 1000g – REF SD400140

This product was developed for use in the dental sector and should be used in accordance with the Instructions for Use. The manufacturer cannot be made liable for any alternative use. In addition, the user is obliged to check the material at their own responsibility for its suitability and application for the intended purposes, especially if these are not listed in the Instructions for Use. Serious incidents must be reported to DeltaMed GmbH as well as the competent authorities.

Only

Important product information

Federal (U.S.) law restricts the finished device to sale by or on the order of a dentist.

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SHERA Werkstoff-Technologie GmbH & Co. KG