voco
V-Print® model 2.0
EN Instructions for use

Product description:
V-Print model 2.0 is a light-curing resin for the generative manufacture of dental models using CAD/CAM technology.

Colour:
Beige

Indication:
Dental working and presentation models

Contraindications:
V-Print model 2.0 contains methacrylate and phosphine oxide. In case of any known hypersensitivities (allergies) to these ingredients of **V-Print model 2.0**, the application should be dispensed with.
Not suitable for intraoral application.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
V-Print model 2.0 should only be used by a professionally trained dental practitioner.

CAD software¹ dental scanner	Software for the planning and design of dental models. The software and dental scanner must satisfy local and current standards and the manufacturer's instructions for the assurance of the patient-specific design as an STL data set.
	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created. For example: - Autodesk Netfabb version 2020 or later for SoliFlex 3D printing.

The designation **Software as Medical Device SMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SoliFlex 170 (PowerVat) VOCO SoliFlex 350 (PowerVat) VOCO SoliFlex 650 (PowerVat) VOCO SoliFlex 170 HD
Post-curing devices	For example: Ofoflash G171

See also: accompanying list of resources or www.voco.dental/3dprintingpartners
All manuals and/or operating instructions for the respective programmes, and for devices, devices and parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programmes, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

CAUTION: Changes to the process equipment, parameters, or software may result in the end object from **V-Print model 2.0** not meeting the specifications.

Use:

Preparation:
For an indication-appropriate CAD construction, the following design conditions must be observed:

- Recommended wall thickness if the model is not filled: 3 mm

Prepare a print job using slicing software, in addition to the material-dependent construction specifications in these instructions for use, please also observe the dependencies of the positioning, support type, and fit found in our other publications, for your construction. The **V-Print model 2.0** has been developed for a high-resolution application. It is thus recommended to select a low layer thickness when generating the print data set.

Processing:

Note: Use separate material containers and cleaning baths for each printing material, in order to avoid cross-contamination.

Note: Shaking the material before printing is not necessary.

The materials container should be filled immediately before the start of the printing process. It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.

Start the print job observing the parameters that you previously selected.

Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, the fully printed objects must be removed from the build platform, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.

Once your work is completed, transfer the remaining material from the materials container into the original container. This allows for the materials container to be inspected and facilitates optimal storage of the printing material.

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